



PHILIPS

Patient Monitor


G30E

**Compact, versatile
and cost-effective**
Philips G30E B66 Option

Supporting your caregivers

A patient's condition can change rapidly. Detecting deterioration and responding to fluctuations in vital signs is a key part of care staff's role. Having a patient monitoring solution that delivers fast reliable insights can help save precious time and support informed clinical decision making.

Designed and built in line with Philips high standards of quality, the G30E B66 Option patient monitor provides valuable insights into essential parameters and is simple for your caregivers to adopt and operate.



Ergonomic design, market proven capabilities

The monitor's lightweight, ergonomic design makes it easy and practical to use. The device includes essential measurements, such as the oxyCRG and Mini-trend software functions that provide basic clinical support. The G30E B66 Option complies with the European Union Medical Device Regulation (MDR).

Each monitor can be connected to the Philips Efficia CMS200 Central Station and includes an HL7 data output function for EMR/HIS integration. Moreover, the solution features the Philips FAST SpO₂ and enhanced Philips ECG (ST/AR) algorithms. These essential market-proven capabilities support accurate measurements, helping you to make informed clinical decisions.

A closer look at the G30E B66 Option

Lightweight, versatile and ready to meet your monitoring needs

Fanless design to reduce dust deposits and control cross-infection risks

With the Philips ECG and FAST SpO₂ algorithm

Rechargeable lithium ion battery

3-grade alarm LED for color-coded visual alarms

Easy-to-read, 10.4-inch display with large, color-coded numerics and waveforms

Navigation wheel for information input and screen navigation

Backlit function buttons for quick access to key functions



G30E B66 Option at a glance



Ergonomic design and easy to move around with a weight of just 3.5 kilos.



Bed-to-bed overview function including waveforms for up to 32 monitors on the same local network



Delivers accurate measurements with an extensive range of market-proven functions and features.



Allows 2.5 hours of continuous monitoring thanks to 4800 mAh, 11.1V rechargeable lithium ion battery



Facilitates a low total cost of ownership with affordable servicing, parts and accessories



HL7 data output function for EMR and HIS integration

Reliable, familiar consumables

Philips G30E B66 Option is designed for use with the same consumables used on the Philips Efficia CM series monitors. This means you can share the products across devices and care settings, helping lower your total cost of ownership.

These durable, reliable consumables have been designed in line with Philips market proven quality standards.



High quality with a low total cost of ownership

The G30E B66 Option monitor delivers proven Philips technology at an attractive price point. This makes it an ideal option for healthcare facilities seeking to deliver outstanding patient care across a range of acuity levels with a manageable level of capital investment and low total cost of ownership.

Discover this easy-to-learn, preconfigured bedside monitor, with proven Philips technology today. For more information, please contact your Philips representative.

To learn more about Philips patient monitor G30E, please reach out to your Philips representative.

